

PVD reflected light microscope

Artikel-Nr. 33229

Modern microscopy made easy!
This reflected light microscope is ideal for determining coating thicknesses based on the calotte grinding method and is therefore the perfect companion to our Calotte L.

Experten kontaktieren

+49 2405 44 70 122

coatingtechnology@cemecon.de



Technische Daten

Type	BX53M
Coating procedure	Pretreatment Follow up treatment
Process step	Incoming inspection Quality control
Test method	Microscopy PVD coating thickness
Coating Technology	PVD SPUTTER PVD HIPIMS

The large magnification allows evaluation of the layer structure and the layer/support material interface. The BX53M simplifies complex microscopy tasks with its well-designed and user-friendly controls. Users can get the most out of the microscope without extensive training. The simple and convenient operation of the BX53M also improves reproducibility by minimizing human error.

Your advantages:

- Simple light source
- Intuitive microscope controls
- Fast focus adjustment
- Consistent lighting
- Easy and ergonomic operation
- Easy restoration of microscope settings
- Basic measurements

Scope of delivery:

- Reflected light microscope
- Included image processing „PRECiV Pro" Version 1.2
- Included PRECiV „Layer Thickness"
- Camera, Workstation, Monitor